

Anti-WDR3 antibody (60-300 aa) (STJ11106489)
STJ11106489

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IF/ICC:1:50-1:200
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	10885
Gene Symbol	WDR3
Uniprot ID	WDR3_HUMAN
Immunogen	
Immunogen Region	60-300 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 60-300 of human WDR3 (NP_006775.1).
Immunogen Sequence	ILQGLKQEVTCLCPSDGLH LAVGYEDGSIRIFSLSGEG NVTFNHGKAAITLKYDQLG GRLASGSKDIDIIVWDVINE SGLYRLKGHKDAITQALFLR EKNLLVTSGKDTMVKWWDLDTQHCFKTMVGHRTVWGLVL LSEEKRLITGASDSELRVWD IAYLQEIETPEEPDPKKIKG SSPGIQDTLEAEDGAFETDE APEDRILSCRKAGSIMREGR DRVVNLAVDKTGRILACHG