

Anti-DDHD2 antibody (360-500 aa) (STJ11106493)

STJ11106493

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	23259
Gene Symbol	DDHD2
Uniprot ID	DDHD2_HUMAN
Immunogen	
Immunogen Region	360-500 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 360-500 of human DDHD2 (NP_056029.2).
Immunogen Sequence	ILTNQKDSLGLDIDSEKDSLNI VMDQGDPTLEEDLKKLQL SEFFDIFEKEKVDKEALALC TDRDLQEIGIPLGPRKKILN YFSTRKNSMGIKRPAPQPAS GANIPKESEFCSSSNTRNGD YLDVGIGQVSVKYPRLIYKP E