

Anti-UBXN7 antibody (200-390 aa) (STJ11106494)

STJ11106494

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	26043
Gene Symbol	UBXN7
Uniprot ID	UBXN7_HUMAN
Immunogen	
Immunogen Region	200-390 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 200-390 of human UBXN7 (NP_056377.1).
Immunogen Sequence	HFIFWQVYHDSEEGQRYIQF YKLGDFPYVSILDPRTGQKL VEWHLQDVSSFLDQVTGFLG EHGQLDGLSSSPKKCARSE SLIDASEDSQLEAAIRASLQ ETHFDSTQTKQDSRSDEESE SELFSGSEEFISVCGSDEEE EVENLAKSRKSPHKDLGHRK EENRRPLTEPPVRTDPGTAT NHQGLPAVDSE