

Anti-UNC45A antibody (1-200 aa) (STJ11106498)

STJ11106498

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	55898
Gene Symbol	UNC45A
Uniprot ID	UN45A_HUMAN
Immunogen	
Immunogen Region	1-200 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human UNC45A (NP_061141.2).
Immunogen Sequence	MTVSGPGTPEPRPATPGASS VEQLRKEGNELFKGDYGGG LAAYTQALGLDTPQDQAVL HRNRAACHLKLEDYDKAETE ASKAIEKDGGDVKALYRRSQ ALEKLGRDQAVLDLQRCVS LEPKNKVFQEALRNIGGQIQ EKVRVMSSTDAKVEQMFQIL LDPEEKGTETKKQKASQNLVV LAREDAEAEKIFRSNGVQLL