

Anti-GARNL3 antibody (250-500 aa) (STJ11106506)
STJ11106506

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	84253
Gene Symbol	GARNL3
Uniprot ID	GARL3_HUMAN
Immunogen	
Immunogen Region	250-500 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 250-500 of human GARNL3 (NP_115669.3).
Immunogen Sequence	RGGLDTKNDTTGIHSVYTVY QGHEIMFHVSTMLPYSKENK QQVERKRHIGNDIVTIVFQE GEESSPAFAKPSMIRSHFTHI FALVRYNQNDNYRLKIFSE ESVPLFGPPLPTPPVFTDHQ EFRDPLLVLINGEKATLET PTFAQKRRTLDMLIRSLHQ DLMPDLHKMNLNRRSFSDVL PESPKSARKKEEARQAEFVR IGQALKLKSIVRGDAPSSLA ASGICKKEPWEPQCFCSNF