

Anti-WDR83 antibody (1-315 aa) (STJ11106509)

STJ11106509

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IF/ICC:1:50-1:200
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	84292
Gene Symbol	WDR83
Uniprot ID	WDR83_HUMAN
Immunogen	
Immunogen Region	1-315 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-315 of human WDR83 (NP_115708.1).
Immunogen Sequence	MAFPEPKPRPELPQKRLKT LDCGQGAVRAVRFNVDGNYC LTCGSDKTLKLVNPLRGTLT RTYSGHGYEVLDAAGSFDNS SLCSGGGDKAVVLWDVASGQ VVRKFRGHAGKVNTVQFNEE ATVILSGSIDSSIRCWDCRS RRPEPVQTLDEARDGVSSVK VSDHEILAGSVDGRVRRYDL RMGQLFSDYVGSPITCTCFS RDGQCTLVSSLDSTLRLLDK DTGELLGEYKGHKNQEYKL