

Anti-ACOT4 antibody (180-340 aa) (STJ11106512)

STJ11106512

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IHC-P:1:50-1:200
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	122970
Gene Symbol	ACOT4
Uniprot ID	ACOT4_HUMAN
Immunogen	
Immunogen Region	180-340 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 180-340 of human ACOT4 (NP_689544.3).
Immunogen Sequence	GHGFATLALAYNFEDLPNN MDNISLEYFEEAVCYMLQHP QVKGPGIGLLGISLGADICL SMASFLKNVSATVSINGSGI SGNTAINYKHSSIPPLGYDL RRIKVAFGSLVDIVDIRNAL VGGYKNPSMIEKAQGPIL LIVGQDDHNWRSELYAQTVS E