

Anti-PUSL1 antibody (150-250 aa) (STJ11106514)

STJ11106514

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	126789
Gene Symbol	PUSL1
Uniprot ID	PUSL1_HUMAN
Immunogen	
Immunogen Region	150-250 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 150-250 of human PUSL1 (NP_699170.1).
Immunogen Sequence	NLCWTLPADCLDMVAMQEAA QHLLGTHDFSAFQSAGSPVP SPVRTLRRVSVSPGQASPLV TPEESRKLRFWNLEFESQSF LYRQVRRMTAVLVAVGLGAL A