

Anti-PPP1R1C antibody (1-109 aa) (STJ11106516)

STJ11106516

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	151242
Gene Symbol	PPP1R1C
Uniprot ID	PPR1C_HUMAN
Immunogen	
Immunogen Region	1-109 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-109 of human PPP1R1C (NP_001074014.1).
Immunogen Sequence	MEPNSPKKIQFAVPVFSQI APEAAEQIRKRRPTASLVI LNEHNPEIDDKRGPNQTGE LQNASPKQRKQSVYTPPTIK GVKHLKGQNESAFPEEEEGT NEREEQRDH