

Anti-ACR antibody (50-150 aa) (STJ11106527)

STJ11106527

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	49
Gene Symbol	ACR
Uniprot ID	ACRO_HUMAN
Immunogen	
Immunogen	50-150 aa
Region	
Specificity	A synthetic peptide corresponding to a sequence within amino acids 50-150 of human ACR (NP_001088.2).
Immunogen	QHGAWPWMVSLQIFTYN SHR YHTCGGSLLSRWVLTAAHC FVGKNNVHDWRLVFGAKEIT YGNNKPVKAPLQERYVEKII
Sequence	IHEKYN SATEGNDIALVEIT P