

Anti-PRDM16 antibody (75-226 aa) (STJ11106528)

STJ11106528

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	63976
Gene Symbol	PRDM16
Uniprot ID	PRD16_HUMAN
Immunogen	
Immunogen Region	75-226 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 75-226 of human PRDM16 (NP_071397.3).
Immunogen Sequence	PEDIPIPADFELRESSIPGA GLGVWAKRKMEAGERLGPCV VVPRAAAKETDFGW EQILTD VEVSPQEGCITKISED LGSE KFCVDANQAGAGSWLKYIRV ACSCDDQNLTMCQISEQIYY KVIKDIEPGEELLVHVKEGV YPLGTVPPGLDE