

Anti-SLC5A8 antibody (461-610 aa) (STJ11106533)

STJ11106533

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	160728
Gene Symbol	SLC5A8
Uniprot ID	SC5A8_HUMAN
Immunogen	
Immunogen Region	461-610 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 461-610 of human SLC5A8 (NP_666018.3).
Immunogen Sequence	QIYPPLPERTLPLHLDIQGC NSTYNETNLMTTTEMPFTTS VFQIYNVQRTPLMDNWYSLS YLYFSTVGTLLVLLVGILVS LSTGGRKQNLDPRIYLTKEF FLSNFDIFKKKKHVLVSYKSH PVEDGGTDNPAFNHIELNSD QSGKSNGTRL