

Anti-Methyl-hnRNP K-H187 antibody (STJ11106539)

STJ11106539

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	3190
Gene Symbol	HNRNPK
Uniprot ID	HNRPK_HUMAN
Immunogen	CPHST
Immunogen Region	
Specificity	A synthetic methylated peptide around H187 of human hnRNP K (NP_112552.1).
Immunogen Sequence	CPHST