

Anti-Methyl-SUPT5H/SPT5-H418 antibody (STJ11106544)

STJ11106544

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	6829
Gene Symbol	SUPT5H
Uniprot ID	SPT5H_HUMAN
Immunogen	REHNF
Immunogen Region	
Specificity	A synthetic methylated peptide around H418 of human SUPT5H/SPT5 (NP_003160.2).
Immunogen Sequence	REHNF