

Anti-DNA-RNA Hybrid antibody [SMM] (STJ11106552)

STJ11106552

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	DB/ELISA
Host/Source	Mouse
Reactivity	Human/Mouse/Rat/Other

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SMM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	DB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG2ak
Storage Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	Synthetic DNA/RNA
Immunogen Sequence	