

Anti-DNA Ligase III/LIG3 antibody (730-1009 aa) (STJ11106570)

STJ11106570

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	
Applications	WB/IP/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IP:0.5 Mu g-4 Mu g antibody for 200 Mu g-400 Mu g extracts of whole cells ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	3980
Gene Symbol	LIG3
Uniprot ID	DNL3_HUMAN
Immunogen	
Immunogen Region	730-1009 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 730-1009 of human DNA Ligase III/LIG3 (NP_039269.2).
Immunogen Sequence	AGGHDDATLARLQNELDMVK ISKDPSKIPSWLVNKIYYP DFIVDPKKAHVWEITGAEF SKSEAHTADGISIRFPRCTR IRDDKDWKATNLPQLKELY QLSKEKADFTVAGDEGSST TGGSSSEENKGPSSGSAVSRKA PSKPSASTKKAEGKLSNSNS KDGNMQTAKP SAMKVGEKLA TKSSPVKVGEKRAAETLC QTKVLLDIFTGVRLYLPST PDFSRLRRYFVAFDGLVQ



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081