

Anti-VCAM1 antibody (640-739 aa) [SRM] (STJ11106572)

STJ11106572

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:2000-1:12000 IHC-P:1:1000-1:4000 IF/ICC:1:100-1:400 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	7412 NA
Gene Symbol	VCAM1 Vcam1
Uniprot ID	VCAM1_HUMAN VCAM1_MOUSE
Immunogen	640-739 aa
Immunogen Region	
Specificity	A synthetic peptide corresponding to a sequence within amino acids 640-739 of human VCAM1 (P19320).
Immunogen Sequence	ETGDTVLKSIDGAYTIRKAQ LK DAGVYECESKNKVGSQLR SLTLDVQGRENNKDYFSP EL LVL YFASSLIIPAIGMIIYF ARKANMKGSYSLVEAQKSKV