

Anti-HEXB antibody (277-556 aa) (STJ11106573)

STJ11106573

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	3074
Gene Symbol	HEXB
Uniprot ID	HEXB_HUMAN
Immunogen	
Immunogen Region	277-556 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 277-556 of human HEXB (NP_000512.1).
Immunogen Sequence	YARLRGIRVLPEFDTPGHTL SWGKGQKDLLTPCYSRQNKL DSFGPINPTLNTTYSFLTTF FKEISEVFPDQFIHLGGDEV EFKCWESNPQIDFMRQKGF GTDFKKLESFYIQKVLIIA TINKGSIVWQEVFDDKAKLA PGTIVEVWKDSAYPEELSRV TASGFPVILSAPWYLDLISY GQDWRKYKVEPLDFGGTQK QKQLFIGGEACLWGEYVDAT NLTPRLWPRASAVGERLWS