

Anti-LRP2 antibody (4050-4350 aa) (STJ11106580)

STJ11106580

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IF/ICC:1:50-1:200
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	4036
Gene Symbol	LRP2
Uniprot ID	LRP2_HUMAN
Immunogen	
Immunogen Region	4050-4350 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 4050-4350 of human LRP2 (NP_004516.2).
Immunogen Sequence	EGSSPLLLLPDNRIRKYNL SSERFSEYLQDEEYIQAVDY DWDPKDIGLSVYYTVRGEG SRFGAIKRAYIPNFESGRNN LVQEVDLKLYVMQPDGIAV DWVGRHIYWSDVKNKRIEVA KLDGRYRKWLSTDLQPA AIAVNPKLGLMFWTDWGKEPK IESAWMNGEDRNILVFEDLG WPTGLSIDYLNNDRIYWSDF KEDVIETIKYDGTDRRVIK EAMNPYSLDIFEDQLYWIS