

Anti-CD298/ATP1B3 antibody (180-279 aa) [SRM] (STJ11106584)

STJ11106584

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:1500-1:3000 IHC-P:1:100-1:500 IF/ICC:1:200-1:800 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	483
Gene Symbol	ATP1B3
Uniprot ID	AT1B3_HUMAN
Immunogen	
Immunogen Region	180-279 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 180-279 of human ATP1B3 (P54709).
Immunogen Sequence	GLKPEGVPRIDCVSKNEDIP NVAVYPHNGMIDLKYPYYG KKLHVGYLQPLVAVQVSFAP NNTGKEVTVECKIDGSANLK SQDDRDKFLGRVMFKITARA