

Anti-SARS-CoV-2-M₁ protein antibody (113-222 aa) (STJ11106601)

STJ11106601

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	SARS-CoV-2

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	113-222 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 113-222 of coronavirus M ₁ protein (YP_009724393.1).
Immunogen Sequence	NPETNILLNVPLHGTILTRP LLESELVIGAVILRGHLRIA GHHLGRCDIKDLPKEITVAT SRTLSSYYKLGASQQRVAGDSG FAAYSRYRIGNYKLNTDHSS SSDNIALLVQ