

Anti-SARS-CoV-2-ORF8 antibody (1-121 aa) (STJ11106602)

STJ11106602

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	SARS-CoV-2

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	1-121 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-121 of coronavirus ORF8 (YP_009724396.1).
Immunogen Sequence	MKFLVFLGIITTVAAAFHQEC SLQSQCTQHQPYYVDDPCPIH FYSKWYIRVGARKSAPLIEL CVDEAGSKSPIQYIDIGNYT VSCLPFTINCQEPKLGSLVV RCSFYEDFLEYHDVRVVLDF I