

Anti-ATP6V1F antibody (1-119 aa) (STJ11106608)

STJ11106608

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	9296
Gene Symbol	ATP6V1F
Uniprot ID	VATF_HUMAN
Immunogen	
Immunogen Region	1-119 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-119 of human ATP6V1F (NP_004222.2).
Immunogen Sequence	MAGRGLIAVIGDEDVTGF LLGGIGELNKNRHPNFLVVE KDTTINEIEDTFRQFLNRDD IGIIINQYIAEMVRHALDA HQQSIPAVLEIPSKEHPYDA AKDSILRRARGMFTAEDLR