

Anti-PIEZO2 antibody (2480-2658 aa) (STJ11106612)

STJ11106612

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	63895
Gene Symbol	PIEZO2
Uniprot ID	PIEZ2_HUMAN
Immunogen	
Immunogen Region	2480-2658 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 2480-2658 of human PIEZO2 (NP_071351.2).
Immunogen Sequence	TGAMQFLENYEKEDITVAEL EGNNSLWTISPPSKQKMIH ELLDPNSSFSVVFWSIQRN LSLGAKSEIATDKLSFPLKN ITRKNIAKMIAGNSTESSKT PVTIEKIYPYYVKAPSDSNS KPIKQLSENFMFDITILS RDNTTKYNSEWWVLNLTGNR IYNPNSQALELVFNDKVS