

Anti-Fam221a antibody (1-292 aa) (STJ11106613)
STJ11106613

GENERAL INFORMATION




Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	231946
Gene Symbol	Fam221a
Uniprot ID	F221A_MOUSE
Immunogen	
Immunogen Region	1-292 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-292 of mouse Fam221a (NP_766315.1).
Immunogen Sequence	MERLTLPGGAEAVDEYLEY RRIVGEDDGGKLFTEEYEE YKKKVLPMRLQNRLFVSWRS PTGMDCKLVGPETLCFCTHR YKQHKTDFTIPQQRPIALP CRVSGCRCKAYHYVPLNGTQ PIRCRCKHFADQHSALGFT CNACSKCSGFHSCYTCACGQ PAYAHDTVVETRQERLAQ GK PVGRDVPYAAMGGLTG FSSL AEGYMRLLDDSGIGAPSV EVL DSAVSAMDHPFLRAMHAPS



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081