

Anti-TYRP1 antibody (350-450 aa) (STJ11106625)

STJ11106625

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	7306
Gene Symbol	TYRP1
Uniprot ID	TYRP1_HUMAN
Immunogen	
Immunogen Region	350-450 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 350-450 of human TYRP1 (NP_000541.1).
Immunogen Sequence	NSTNSFRNTVEGYSDPTGKY DPAVRS LHNLAHLFLNGTGG QTHLSPNDPIFVLLHTFTDA VFDEWLRRYNADISTFPLEN APIGHN RQYNMVPFWPPVTN T