

Anti-ASH1L antibody (2074-2293 aa) (STJ11106636)
STJ11106636

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:1000-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	55870
Gene Symbol	ASH1L
Uniprot ID	ASH1L_HUMAN
Immunogen	
Immunogen Region	2074-2293 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 2074-2293 of human ASH1L (NP_060959.2).
Immunogen Sequence	SNVYVDVKPLSGYEATTCNC KKPDDDTTRKGCVDDCLNRMIFAECSPNTCPGGEQCCNQRIQRHEWVQCCLERFRAEEKGWG IRTKPEPLKAGQFIIEYLGEV VSEQEFRNRMIEQYHNHSDH YCLNLDSGMVIDSYRMGNEA RFINHSCDPNCEMQKWSVNG VYRIGLYALKDMPAGTELTY DYNFHSFNVEKQQLCKCGFE KCRGIIGGKSQRVNGLTSSK