

**Anti-alt-RPL36 antibody (1-148AA aa) (STJ11106637)**

STJ11106637

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	
<b>Applications</b>	WB/ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	Lot specific
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purification
<b>Dilution</b>	WB:1:500-1:2000
<b>Range</b>	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
<b>Formulation</b>	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
<b>Isotype</b>	IgG
<b>Storage</b>	
<b>Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	
<b>Immunogen</b>	1-148AA aa
<b>Region</b>	
<b>Specificity</b>	Recombinant fusion protein containing a sequence corresponding to amino acids 1-148AA of alt-RPL36.
<b>Immunogen</b>	MASEAEGAGTRGSGPRELSR YSPPLLPAVVRSHGPTLPY GRGPQQGPQSDQEREQAQAQ PTPRASDQTHQVRAGHDSGG
<b>Sequence</b>	VWLCPVRAARHGVTEGLQGQ TGPQIQEKGGDAHPRQEEA GGAEQRTGRHEESCCQERLS PSPALSLK