

Anti-PIR antibody (1-290 aa) (STJ11106645)

STJ11106645

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	8544
Gene Symbol	PIR
Uniprot ID	PIR_HUMAN
Immunogen	
Immunogen Region	1-290 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-290 of human PIR (NP_003653.1).
Immunogen Sequence	MGSSKKVTLVLSREQSEGV GARVRRSIGRPELKNLDPFL LFDEFKGGRRPGGFPDHPHRG FETVSYLLEGGSM AHEDFCG HTGKMNPGLQWMTAGRGIL HAEMPCSEEP AHGLQLWVNL RSSEKMVEPQYQELKSEEIP KPSKDGVTVAVISGEALGIK SKVYTRTP TLYLDFKLDPGA KHSQPIPKGWTSFIYTISGD VYIGPDDAQQKIEPHHTAVL GEGDSVQVENKDKRSHFV