

Anti-GLRX5 antibody (32-157 aa) (STJ11106647)

STJ11106647

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:1000-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	51218
Gene Symbol	GLRX5
Uniprot ID	GLRX5_HUMAN
Immunogen	
Immunogen Region	32-157 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 32-157 of human GLRX5 (NP_057501.2).
Immunogen Sequence	AGSGAGGGGSAEQLDALVKK DKVVVFLKGTPEQPQCGFSN AVVQILRLHGVRDYAAYNVL DDPQLRQGIKDYSNWPTIPQ VYLNGEFVGGCDILLQMHQN GDLVEELKKLGIHSALLDEK KDQDSK