

Anti-TMPRSS13 antibody (368-567 aa) (STJ11106648)

STJ11106648

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	84000
Gene Symbol	TMPRSS13
Uniprot ID	TMPSD_HUMAN
Immunogen	
Immunogen Region	368-567 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 368-567 of human TMPRSS13 (NP_001070731.1).
Immunogen Sequence	FFVTREKVLGKWKVYAGTSN LHLPEAAASIAEIIINSNYT DEEDDYDIALMRLSKPLTLS AHIIHPACLPMHGQTFSLNET CWITGFGKTRETDDKTSPLF REVQVNLIDFKKCNLYLVYD SYLTPRMMCAGDLRGGRDSC QGDGGGPLVCEQNNRWYLAG VTSWGTGCGQRNKPQVYTKV TEVLPWIYSKMESEVRFKRS