

Anti-RGL2 antibody (1-240 aa) (STJ11106667)

STJ11106667

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	5863
Gene Symbol	RGL2
Uniprot ID	RGL2_HUMAN
Immunogen	
Immunogen Region	1-240 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human RGL2 (NP_004752.1).
Immunogen Sequence	MLPRPLRLLLDTSPPGGVVL SSFRSRDPEEGGPGGLVVG GGQEEEEEEEEAPVSVWDE EEDGAVFTVTSRQYRPLDPL VPMPPPRSSRRRLRAGTLEAL VRHLLDTRTSGTDVSFMSAF LATHRAFTSTPALLGLMADR LEAESHPTDELERTTEVAI SVLSTWLASHPEDFGSEAKG QLDRLSFLLQTGYAAGKGV GGGADLIRNLRSRVDPQAP DLPKPLALPGDPPADPTDV