

Anti-B3GAT3 antibody (29-190 aa) (STJ11106707)

STJ11106707

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	26229
Gene Symbol	B3GAT3
Uniprot ID	B3GA3_HUMAN
Immunogen	
Immunogen Region	29-190 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 29-190 of human B3GAT3 (NP_036332.2).
Immunogen Sequence	QPCDCLPPLRAAAEQLRQKD LRISQLQAELERRPPAPAQP PEPEALPTIYVVTPTYARLV QKAEVLRLSQTLSLVPRLHW LLVEDAEGPTPLVSGLLAAS GLLFTHLWLTTPKAQLREG EPGWVHPRGVEQRNKALDWL RGRGGAVGGEKDPPTPGTQG VV