

Anti-MRPS15 antibody (58-180 aa) (STJ11106715)

STJ11106715

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	64960
Gene Symbol	MRPS15
Uniprot ID	RT15_HUMAN
Immunogen	
Immunogen Region	58-180 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 58-180 of human MRPS15 (NP_112570.2).
Immunogen Sequence	PAQSRLLLLDDPPSTLLKDYQ NVPGIEKVDVVKRLLSLEM ANKKEMLKIKQEQFMKKIVA NPEDTRSLEARIIALSVKIR SYEEHLEKHKDKAHKRYLL MSIDQRKKMLNLRNTNYDV FEK