

Anti-DLG5 antibody (1436-1701 aa) (STJ11106727)

STJ11106727

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	9231
Gene Symbol	DLG5
Uniprot ID	DLG5_HUMAN
Immunogen	
Immunogen Region	1436-1701 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1436-1701 of human DLG5 (NP_004738.3).
Immunogen Sequence	SHSRSSSHLDPAAGTHSTLQG SGTTPPEHPSVIDPLMEQDE GPSTPPAKQSSSRIAGDANK KTLPRVVFIIKKSQLELGVH LCGGNLHGVFVAEVEDDSPA KGPDGLVPGDLILEYGS LDV RNKTVEEVYVEMLKPRDGVR LKVQYRPEEFTKAKGLPGDS FYIRALYDRLADVEQELESFK KDDILYVDDTLPGGTFGSWM AWQLDENAQKIQRGQIPSKY VMDQEFSSRLSMSEVKDDN