

Anti-BTNL9 antibody (35-250 aa) (STJ11106730)

STJ11106730

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	153579
Gene Symbol	BTNL9
Uniprot ID	BTNL9_HUMAN
Immunogen	
Immunogen Region	35-250 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 35-250 of human BTNL9 (NP_001295174.1).
Immunogen Sequence	SSEVKVLGPYPILALVGEE VEFPCHLWPQLDAQQMEIRW FRSQTFNVVHLYQEQQELPG RQMPAFRNRTKLVKDDIAYG SVVLQLHSIIPSDKGTYGCR FHSDNFSGEALWELEVAGLG SDPHLSLEGFKEGGIQLRLR SSGWYPKPKVQWRDHQGQCCL PPEFEAIWDAQDLFSLETS VVVRAGALSNVSVSIQNLLL SQKKELVWQIADVFP