

Anti-ADGRG3 antibody (21-270 aa) (STJ11106736)

STJ11106736

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	222487
Gene Symbol	ADGRG3
Uniprot ID	AGRG3_HUMAN
Immunogen	
Immunogen Region	21-270 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 21-270 of human ADGRG3 (NP_740746.4).
Immunogen Sequence	QEKPTGPRNTCLGSNNMYD IFNLNDKALCFTKCRQSGSD SCNVENLQRYWLNYYEHLMK EGLTQKVNTPFLLKALVQNLS TNTAEDFYFSLEPSQVPRQV MKDEDKPPDRVRLPKSLFRS LPGNRSVVRLAVTILDIGPG TLFKGPRLGLDGSGLVNNR LVGLSVGQMHVTKLAEPLEI VFHQRPPTMTLTCVFDV TKGTTGDSSEGCSTEV RPE GTVCCDHLTFALLLRPT