

Anti-GPR98 antibody (30-300 aa) (STJ11106737)

STJ11106737

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB: 1:500-1:2000
Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	84059
Gene Symbol	ADGRV1
Uniprot ID	AGRV1_HUMAN
Immunogen	
Immunogen Region	30-300 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 30-300 of human GPR98 (NP_115495.3).
Immunogen Sequence	ETEIRFTGQTEFVNETSTT VIRLIIRIGEPANVTAIVS LYGEDAGDFFDTYAAAFIPA GETNRTVYIACDDDLPEPD ETFIFHLTLQKPSANVKLGW PRTVTVTILSNDNAFGIISF NMLPSIAVSEPKGRNESMPL TLIREKGTYGMVMVTFEVEG GPNPPDEDLSPVKGNIFFPP GRATVIYNLTVLDDDEVPEND EIFLIQLKSVEGGAEINTSR NSIEIIKKNDSPVRFLQS