

Anti-WNT11 antibody (225-354 aa) (STJ11106741)

STJ11106741

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	7481
Gene Symbol	WNT11
Uniprot ID	WNT11_HUMAN
Immunogen	
Immunogen Region	225-354 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 225-354 of human WNT11 (NP_004617.2).
Immunogen Sequence	KGLQELQDVAADLKTRYLSA TKVVHRPMGTRKHLVPKDLD IRPVKDSELVYLQSSPDFCM KNEKVGSHGTQDRQCNKTSN GSDSCDLMCCGRGYNPYTDR VVERCHCKYHWCCYVTCRRC ERTVERYVCK