

Anti-USH2A antibody (390-520 aa) (STJ11106757)

STJ11106757

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000
	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	7399
Gene Symbol	USH2A
Uniprot ID	USH2A_HUMAN
Immunogen	
Immunogen Region	390-520 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 390-520 of human USH2A (NP_996816.2).
Immunogen Sequence	FSQPQTEIRIQRKENS LDW EDWQYFARNCGAFGMKNNGD LEKPDSVNCQLQSNFTPYSR GNVTF SILTPGPNYRPGYNN FYNTPSLQEFVKATQIRFHF HGQYYTTETAVNLRHRYAV DEITISGRQCQC