

Anti-IL19 antibody (1-215 aa) (STJ11106763)

STJ11106763

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	29949
Gene Symbol	IL19
Uniprot ID	IL19_HUMAN
Immunogen	
Immunogen Region	1-215 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-215 of human IL19 (NP_715639.1).
Immunogen Sequence	MCTEGAFPHRSACSLPLTHV HTHIHVCVPVLWGSVPRGMK LQCVSLWLLGTILILCSV DN HGLRRCLISTDMHHIEESFQ EIKRAIQAKDTFPNVTILST LETLQIIKPLDVCCVTKNLL AFYVDRVFKDHQEPNPKILR KISSIANSFLYMQKTLRQCQ EQRQCHCRQEATNATRVIH D NYDQLEVHAAIKSLGELDV FLAWINKNHEVMFSA