

Anti-SEZ6L2 antibody (600-740 aa) (STJ11106764)

STJ11106764

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	26470
Gene Symbol	SEZ6L2
Uniprot ID	SE6L2_HUMAN
Immunogen	
Immunogen Region	600-740 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 600-740 of human SEZ6L2 (NP_001107572).
Immunogen Sequence	GHRTASDAGFPVGSVQYRC LPGYSLEGAAMLTCSYRDTG TPKWSDRVPKCALKEYEPCLN PGVPENGYQTLYKHYYQAGE SLRFFCYEGFELIGEVTITC VPGHPSQWTSQPPLCKVAYE ELLDNRKLEVTQTDDPSRQL E