

Anti-Prostatic Acid Phosphatase antibody (33-386 aa) [SRM] (STJ11106780)

STJ11106780

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:1000-1:2000
Range	IHC-P:1:1000-1:5000
	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	55
Gene Symbol	ACP3
Uniprot ID	PPAP_HUMAN
Immunogen	
Immunogen Region	33-386 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 33-386 of human Prostatic Acid Phosphatase (ACPP) (NP_001090.2).
Immunogen Sequence	KELKFVTLVFRHGDRSPIDT FPTDPIKESSWPQGFGQLTQ LGMEQHYELGEYIRKRYRKFLNESYKHEQVYIRSTDVDRT LMSAMTNLAALFPPEGVSIW NPILLWQPIPVHTVPLSEDQ LLYLPFRNCPRFQELESETL KSEEFQKRLHPYKDFIATLG KLSGLHGQDLFGIWSKVYDP LYCESVHNFTLPSWATEDTM TKLRELSELSLLSLYGIHKQ KEKSRLQGGVLVNEILNHM