

Anti-human IgG-Fc antibody [SMM] (STJ11106782)
STJ11106782

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IHC-P/ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Clone ID	SMM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC-P:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG1k
Storage Instruction	

TARGET INFORMATION


Gene ID	NA NA NA NA NA NA NA NA NA NA NA NA NA
Gene Symbol	IGHG2 IGHG3 IGHG2_HUMAN IGHG3_HUMAN IGHG4_HUMAN
Immunogen	
Immunogen Region	
Specificity	A synthesized peptide derived from human human IgG (Fc).
Immunogen Sequence	



Immunohistochemistry analysis of paraffin-embedded Human Follicular lymphoma using human IgG (Fc) Mouse mAb (STJ11106782) at dilution of 1:200 (40x lens). High pressure antigen retrieval performed with 0.01M Tris/EDTA Buffer (pH 9.0) prior to immunohistochemistry staining.



Immunohistochemistry analysis of paraffin-embedded Human mantle cell lymphoma using human IgG (Fc) Mouse mAb (STJ11106782) at dilution of 1:200 (40x lens). High pressure antigen retrieval performed with 0.01M Tris/EDTA Buffer (pH 9.0) prior to immunohistochemistry staining.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using human IgG (Fc) Mouse mAb (STJ11106782) at dilution of 1:200 (40x lens). High pressure antigen retrieval performed with 0.01M Tris/EDTA Buffer (pH 9.0) prior to immunohistochemistry staining.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081