

Anti-TPL antibody (200-300 aa) (STJ11106785)

STJ11106785

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	A.thaliana

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	200-300 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 200-300 of arabidopsis thaliana TPL (NP_849672.1).
Immunogen Sequence	VDHSCGPPNGARAPSPVNNP LLGGIPKAGGFPPLGAHGPF QPTASPVPTPLAGWMSSPSS VPHPAVSAGAIALGGPSIPA ALKHPRTPPTNASLDYPSAD S