

Anti-PHYA antibody (836-1122 aa) (STJ11106787)

STJ11106787

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	A.thaliana

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	836-1122 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 836-1122 of arabidopsis thaliana PHYA (NP_172428.1).
Immunogen Sequence	KVSFAFFTRGGKYVECLLCV SKKLDREGVVTGVFCFLQLA SHELQQALHVQRLAERTAVK RLKALAYIKRQIRNPLSGIM FTRKMIEGTELGPEQRRILQ TSALCQKQLSKILDDSDLES IIEGCLDLEMKEFTLNEVLT ASTSQVMMKSNKSVRITNE TGEVMSDTLYGDSIRLQQV LADFMLMAVNFTPSGGQLTV SASLRKQDLGRSVHLANLEI RLTHTGAGIPEFLLNQMFQ