

Anti-GL2 antibody (492-746 aa) (STJ11106789)
STJ11106789

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	A.thaliana

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	492-746 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 492-746 of arabidopsis thaliana GL2 (NP_565223.1).
Immunogen Sequence	ATNVPTKDSLGVTTLAGRKS VLKMAQRMTQSFYRAIAASS YHQWTKITTKTGQDMRVSSR KNLHDPGEPTGVIVCASSSL WLPVSPALLDFFRDEARRH EWDALSNGAHVQSIANLSKG QDRGNSVAIQTVKSREKSIW VLQDSSSTNSYESVVVYAPVD INTTLVLVAGHDPSNIQILP SGFSIIPDGVESRPLVITST QDDRNSQGGSLTLALQTLI NPSPAAKLNMESVESVTNL