

Anti-ATG4A antibody (289-398 aa) (STJ11106792)

STJ11106792

GENERAL INFORMATION




Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	115201
Gene Symbol	ATG4A
Uniprot ID	ATG4A_HUMAN
Immunogen	
Immunogen Region	289-398 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 289-398 of human ATG4A (NP_443168.2).
Immunogen Sequence	TEENGTVNDQTFHCLQSPQR MNILNLDPSVALGFFCKEEK DFDNWCSLVQKEILKENLRM FELVQKHPSHWPPFVPPAKP EVTTTGAEFIDSTEQLEEFD LEEDFEILSV



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081