

Anti-Rad52 antibody (125-418 aa) (STJ11106807)

STJ11106807

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB: 1:500-1:2000
Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	5893
Gene Symbol	RAD52
Uniprot ID	RAD52_HUMAN
Immunogen	
Immunogen Region	125-418 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 125-418 of human Rad52 (NP_602296.2).
Immunogen Sequence	GYGVSEGLKSKALSLEKARK EAVTDGLKRALRSFGNALGN CILDKDYLRSLNKLPRQLPL EVDLTAKARQDLEPSVEEAR YNSCRPNMALGHPQLQQVTS PSRPSHAVIPADQDCSSRSL SSSAVESEATHQRKLRQKQL QQQFRERMEKQQVRVSTPSA EKSEAAPPAPVTHSTPVTV SEPLLEKDFLAGVTQELIKT LEDNSEKWAVTPDAGDGVVK PSSRADPAQTSDTLALNNQ