

Anti-MTMR14 antibody (1-280 aa) (STJ11106808)

STJ11106808

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	64419
Gene Symbol	MTMR14
Uniprot ID	MTMRE_HUMAN
Immunogen	
Immunogen Region	1-280 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human MTMR14 (NP_001070993.1).
Immunogen Sequence	MAGARAAAAAASAGSSASSG NQPPQELGLGELLEFSRTQ YRAKDGGSGTGGSKVERIEKR CLELFGRDYCFSVIPNTNGD ICGHPYRHHVLEYESSEKE KDTFESTVQVSKLQDLIHS KMARCRGRFVCPVILFKGKH ICRSATLAGWGELYGRSGYN YFSGGADDAWADVDEVTTEE DCALRSGDTHLFDKVRGYDI KLLRYLSVKYICDLMVENKK VKFGMNVTSSEKVDKAQRY



Western blot analysis of various lysates using MTMR14 Rabbit pAb (STJ11106808) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution.
Lysates/ proteins: 25 Mu g per lane.
Blocking buffer: 3 % nonfat dry milk in TBST.
Detection: ECL Basic Kit.
Exposure time: 90s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081